

# Annotation Tool

# Starting the Tool

In command prompt:

```
python main.py <group_id>
```

# Two Modes for Contouring

- Nudge mode:  
To make slight changes to the already predicted contour
- Draw Mode:  
To draw new contour using drawing tool

INNER  
Draw mode

Inner/outer wall  
of myocardium

Draw mode/  
Nudge mode

Contour


Navigation  
Instructions

'Q/A' – increase/reduce circle size  
'D' – enter/exit drawing mode  
'S' – save contour  
'N/P' – next/previous image  
'Z/X' – zoom in/out  
'Esc' – quit

MB18044509\_1/series\_CINE\_segmented\_SAX\_b3\_6/image\_1\_inx.txt

# Navigation Instructions

1. 'Q/A' - increase/reduce the size of nudging circle



# Navigation Instructions

2. 'D' - enter/exit drawing mode

Works as line drawing tool in Ms Paint

3. 'S' - save contour
4. 'N/P' - next/previous image
5. 'Z/X' - zoom in/out
6. 'Esc' - quit application

Asks confirmation if contour

is not saved while moving on

